

## Solubilization of membrane proteins using the synthetic polymer diisobutylene-maleic acid (DIBMA)

### Procedure

1. With DIBMA you can solubilize membrane proteins in
  - a. Suspensions of whole cells
  - b. supernatant 9,000
  - c. Homogenized 100,000 Pellet
2. You can add directly the isolated membrane protein solution (40 mg/ml (w/v)) (NOTE: pH is already adjusted with TRIS or HEPES) to your DIBMA/Buffer powder to a final concentration of 0,4 % - 4 %. We advise you to add 2 ml of your solution to 50 mg lyophilized DIBMA to start with a DIBMA concentration of 2,5 %. For screening purposes you can try different DIBMA concentrations to find optimized conditions for solubilization.

DIBMA	4 %	3,5 %	3 %	2,5 %	2 %	1,5 %	1 %	0,8 %	0,6 %	0,4 %
Protein solution	1,25 ml	1,4 ml	1,7 ml	2 ml	2,5 ml	3,33 ml	5 ml	6,25 ml	8,33 ml	12,5 ml

4. Optional: It is possible to screen the total protein concentration (beginning from 40 mg/ml) against 2,5 % DIBMA. The DIBMA or protein concentration screening can influence the success of the solubilization of your target protein.
5. With gentle shaking overnight at 4 °C and/or room temperature you can solubilize your proteins. The optimal conditions are protein specific (screen!).
6. The contained TRIS/HEPES buffer keeps the pH at 7.5 which stabilize the functionality of DIBMA and is suitable for the most applications. If you need a different pH or buffer for your protein solubilization we provide DIBMA without buffer. Please note that a pH smaller than 6.5 is not suitable for solubilization with DIBMA.

## SDS-PAGE Protocol:

If you want to run a SDS Page your proteins have to be separated from the synthetic polymers. The presence of the polymers can lead to smearing bands on SDS-PAGE gels (1,2). Wessel and Flügge discribed the protocol in 1984 originally (3). It was modified in 2017 for the use with DIBMA (1).

1. Measure volume of solubilized protein - polymer sample
2. Vortex the sample with 4x volume of cold methanol
3. Vortex with 1x volume of chloroform
4. Centrifuge for 3 min at 4 °C and 15,000 g with 3x volume of cold water
5. Add 4x volume of methanol to the organic layer and discard the aqueous.
6. Centrifuge at 5,000 g for 1 min to pellet the proteins.
7. Centrifuge at 20,000 g for 5 min at 4 °C.

### References:

1. Oluwole, Abraham Olusegun, et al. „Solubilization of Membrane Proteins into Functional Lipid Bilayer Nanodiscs Using a Diisobutylene/Maleic Acid Copolymer.“ *Angewandte Chemie International Edition* 56.7 (2017): 1919-1924.
2. Lee, Sarah C., et al. „A method for detergent-free isolation of membrane proteins in their local lipid environment.“ *Nature protocols* 11.7 (2016): 1149.
3. Wessel, D. M., and U. I. Flügge. „A method for the quantitative recovery of protein in dilute solution in the presence of detergents and lipids.“ *Analytical biochemistry* 138.1 (1984): 141-143.

